

Prompt Gamma Emission Prediction using an LSTM Network

ECMP 2024 | Munich | 13.09.2024

Fan Xiao¹, Michael Kriechbaum², Domagoj Radonic¹, Niklas Wahl^{3,4}, Ahmad Neishabouri^{4,5}, Nikolaos Delopoulos¹, Katia Parodi², Stefanie Corradini¹, Claus Belka^{1,7,8}, Christopher Kurz¹, Guillaume Landry¹, George Dedes²

¹Department of Radiation Oncology, LMU University Hospital, LMU Munich, Munich, Germany

²Department of Medical Physics, LMU Munich, Munich, Germany

³Department of Medical Physics in Radiation Oncology, German Cancer Research Center (DKFZ), Heidelberg, Germany

⁴National Center for Radiation Oncology (NCRO), Heidelberg Institute for Radiation Oncology (HIRO), Heidelberg, Germany.

⁵Clinical Cooperation Unit Radiation Oncology, German Cancer Research Center (DKFZ), Heidelberg, Germany

⁶German Cancer Consortium (DKTK), partner site Munich, a partnership between DKFZ and LMU University Hospital, Munich, Germany


⁷Bavarian Cancer Research Center (BZKF), Munich, Germany

Background

Proton therapy & Prompt gamma(PG) imaging

■ Proton therapy



- Superior dose conformity compared with photon therapy
 - Limited dose range/reduced exposure to normal tissues
- Limitations due to **range uncertainty**
 - Non-optimal dose coverage of targets/overdose of healthy tissues



Influence of a density heterogeneity to depth dose curves^[1]

■ Prompt gamma (PG) imaging

- PG signals are highly correlated to proton dose distributions
- Potential for **real-time** range monitoring



Schematic setup of PG based proton range verification^[2]


[1] Knopf AC, Lomax A. In vivo proton range verification: a review. *Phys Med Biol.* 2013;58(15):R131-R160.

[2] Tian, Liheng (2020): A new treatment planning concept accounting for prompt gamma imaging for proton range verification. Dissertation, LMU München

Background & Motivation

Fast calculation of expected PG distributions

- The PG fall-off does not explicitly match the therapeutic dose fall-off^[2]




Dose delivery

PG emission

▪ PG emission calculation methods

- MC simulation (slow)
- Analytical dose filtering methods

Fast and accurate prompt gamma calculation



Long short-term memory (LSTM) network-based proton dose calculation^[3]

[2] Tian, Liheng (2020): A new treatment planning concept accounting for prompt gamma imaging for proton range verification. Dissertation, LMU München


[3] Neishabouri A, et al. Long short-term memory networks for proton dose calculation in highly heterogeneous tissues. *Med Phys.* 2021;48(4):1893-1908.

Materials & methods

Patient datasets & MC simulation

■ Patient datasets


- 33 prostate patients treated at LMU University Hospital
- Planning CT scans (resampled voxel size: $1.5 \times 1.5 \times 1.5 \text{ mm}^3$)
- Dataset split: training/validation/testing=20/3/10




■ Geant4 MC simulation for PG and dose distributions

- Convert CT to relative stopping power (RSP) with respect to water
- Number of histories for each PB: **10M**^[4] (20-30 hours)
- Recorded PG threshold: **above 3 MeV**^[5]
- Angular region: $[40^\circ, 140^\circ] \cup [220^\circ, 320^\circ]$, $\Delta\alpha=16.6^\circ$
- PB number simulated per patient: 96
- Cuboid extraction: from patient surface
- Cuboid dimension/resolution: $320 \times 24 \times 24 / 1 \times 2 \times 2 \text{ mm}^3$

MC PG on RSP




MC Dose on RSP




[4] Souris K, et al (2016). Fast multipurpose Monte Carlo simulation for proton therapy using multi-and many-core CPU architectures, *Med Phys*.

[5] Sterpin E, et al (2015). Analytical computation of prompt gamma ray emission and detection for proton range verification, *Phys Med Biol*.



RP (RSP to PG) model: input RSP, output PG



RDP (RSP and dose to PG) model: input RSP and dose, output PG

Materials & methods

Experiments & Metrics

Experiments

- SE: single energy; ME: multi-energy

Cuboid Dataset [energy (PB number)]

Model	Training	Validation	Test
RP	150 MeV (1920)	150 MeV (288)	150 MeV (540)
RDP (SE)	150 MeV (1920)	150 MeV (288)	150 MeV (540)
RDP (ME)	150 MeV (1920), 175 MeV (1920), 200 MeV (1920)	125-200 MeV (288×3)	150 MeV (540), 175 MeV (540), 200 MeV (540), 125-210 MeV (540)

Metrics


- Gamma passing rate γ_{pr} (2%/2mm, PG count > 10% of the maximum PG count)
- Laterally integrated PG profiles with relative count differences ϵ_{rel} (%)
- Range shift $\Delta z^{[6]}$ of two laterally integrated PG-depth profiles

[6] Knopf, A., et al. (2008). Quantitative assessment of the physical potential of proton beam range verification with PET/CT, *Phys Med Biol.* 53(15): 4137.


Results

RP, RDP (SE) and RDP (ME) model comparison


- The RDP(SE) model has smaller Δz than the RP model in the same test dataset (150 MeV)
- Mean γ_{pr} of RDP model (99.1%) is higher than RP model (98.1%)




(a) RP



(b) RDP(SE)



(c) RDP(ME)




Results


RDP (ME) model on four test datasets

- The Δz of test cases from RDP(ME) model is all within 1 mm (except for three outliers)



(a) 129 MeV



(b) 177 MeV



(c) 206 MeV



Results

Runtime

- All computation measurements were carried out on a workstation equipped with an Intel(R) Xeon(R) Gold 6354 3.00 GHz CPU and an NVIDIA RTX A6000 GPU.

Model	Cuboid extraction time (ms)	Model inference time (ms)
RP	55	10
RDP	55x2	11

Conclusions

- A sub-second LSTM-based PG emission prediction method was developed
- Combining RSP and **dose** improves the accuracy of PG range predictions and helps the model generalize to unseen energies

Outlook

- Extend to more treatment sites and spot sizes
- Validation with PG detection data from experimental PG camera needs to be explored



KLINIKUM

Gemeinsam. Fürsorglich. Wegweisend.

Thank you for your attention!



Contact: Fan.Xiao@med.uni-muenchen.de